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IT

FR

POLIVEST is a carbon-free phosphate-bonded precision casting investment for crown and bridgework using all dental alloys by quick or progressive heating procedures.

PHYSICAL DATA (100% liquid concentration)
Mixing ratio (powder/liquid) 100g /24ml
Flow 13 cm
Processing time (22 - 24°C) 5 - 7 min
Setting time (22 - 24°C) 9 - 12 min
Total expansion (linear) 3,00%
Compressive strength 4Mpa

POWDER /LIQUID RATIO
Table with 3 columns: Ring size, Powder, Liquid/dist. water

EXPANSION
Polivest powder is mixed with Polivest mixing liquid. When this liquid is diluted with distilled water, the expansion of the Polivest can be controlled so that compensation can be made for the contraction of the particular casting alloy used.

Table with 5 columns: Type of alloy, Liquid/water ratio, Ring size (x1, x3, x6, x9) for various alloy types like Crown and Bridges, Ceramic Alloy, etc.

STORAGE
Store the powder and liquid at a normal room temperature (23°C). If the liquid is exposed to temperatures at or below 5°C, the liquid will freeze and cannot be used thereafter.

PACKAGES
Polivest powder – 6 kg (40 x 150 g)
Polivest liquid – 1000 ml bottle

DIRECTIONS FOR USE

PREPARATION - Wetting agents are not necessary, but if used ensure that the wax is completely dry before investing. Metal ring method: Use a moistened casting liner.

MIXING - Pre-mix the powder and liquid thoroughly by hand using a spatula. Mix for 60 seconds under a vacuum.

INVESTING - Invest within 6 minutes of the start of mixing (at 23°C). Invest under low vibration. When the ring is full, cease vibration and do not touch the investment until set. Leave to set for 20 minutes from the start of mixing. Scrape the top surface of the investment ring with a sharp knife. Ringless method: After the initial setting (15 min. at 23°C), remove the plastic and allow to set fully.

HEATING PROCEDURES
Quick heating method in a preheated furnace at the final temperature. Moulds (with a metal mould ring) of sizes x1 to x6 can be heated up rapidly. Heat up the furnace according to the type of alloy to be used: 700 – 750°C for gold alloys, 800 – 850°C for ceramic alloys, 900°C for non precious alloys. After 20 minutes setting, place the investment directly into the furnace at the end-temperature.

Heating method in a preheated furnace at 600°C: We recommended that moulds without ring are placed in a preheated furnace at 600°C, just 20 minutes after mixing. After that, heat directly to the final temperature (700°C – 900°C). Hold at the end-temperature for: x 1...40 min x 3...50 min x 6...60 min x 9...90 min

Step heating method: Table with 5 columns: Heating rate, x1, x3, x6, x9

CASTING
Cast in the usual manner: centrifugal, vacuum-pressure casting, etc... Cast as soon as possible after removing the ring from the burnt furnace. Cast according to the processing instructions from the alloy manufacturer.

COOLING
Place the ring upside down. Allow to cool to room temperature.

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