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Trdota po Vickersu: > 26 HV
Upogibna trdnost: > 114 MPa
E-Modul: > 2771 MPa
Prosta monomera: < 1%

DIMENZIJE

Premer: 95, 98, 98.5, 100 mm
Debelina:

- Basic: 12, 13, 14, 15, 16, 18, 20, 25, 30, 35 mm
- Multilayer: 15, 16, 18, 20 mm

POLIDENT PMMA CAD/CAM DISK BASIC (ENOBARVNI) IN MULTILAYER (VEČSLOJNI)

POLIDENT PMMA CAD /CAM DISK tip BASIC (enobarvni) in MULTILAYER so primerni za izdelavo začasnih kron in mostičkov ter za preverjanje naleganja mostičkov na mavčni model ali v ustih pred končnim rezkanjem cirkona.

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POLIDENT PMMA CAD/CAM DISC BASIC (MONOCHROME) AND MULTILAYER

POLIDENT PMMA CAD/CAM DISC BASIC (monochrome) and MULTILAYER are suitable for manufacturing of long-term temporary crowns and bridges and for verifying the bridge on the plaster model/mouth before the final ZrO₂ milling process.

POLIDENT PMMA CAD/CAM DISC MULTILAYER is composed by several different layers. Layers colours are adapted to the natural colour gradient of the tooth. The final result is a very aesthetic and economic restoration.

SOME RULES TO CONSIDER THROUGH THE PMMA C & B DESIGN:

- The minimum wall thickness has to be 1,0mm in occlusal region and 0,6mm in cervical region.
- The connector profile has to be more than 10mm² in anterior region and more than 15 mm² in posterior region.
- A max of two pontics is allowed between two abutment teeth in posterior region.
- Pay attention to the positioning of the restoration within the MULTILAYER disc. The enamel layer is marked on the top of disc.

FINISHING C & B AND CEMENTING IN THE MOUTH:

- Use burs for acrylic. Polishing is possible with polishing paste and soft hair brushes. Pay attention not to overheat the material!
- For cementing the bridge use temporary cement.

SHADES

- Basic: A1 / A2 / A3 / B1 / B2 / B3 / BL1 / BL2 / BL3 / E1 / E2, Clear / Berry / Blue / Green / Orange / DE
- Multilayer: A1 / A2 / A3 / B1

STORAGE

Store at room temperature.

WARNING

In case of known allergies or hypersensitivity to acrylic materials do not use the product. Due to lack of information regarding the safety of PMMA exposure for pregnant women, we advise to avoid using the device for pregnant women.

In case that any serious incident has occurred in relation to the device it should be reported to the manufacturer (Polident d.o.o.; vigilanca@polident.si) and the competent authority of the Member State, where the user is established.

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