GENERAL AND SAFETY INSTRUCTIONS!

- There is a danger of flame occurrence during rapid heating, therefore keep the furnace door closed for at least 15 minutes after the insertion of the cylinder. During this period, do not insert any other cylinders into the furnace.
- Pay attention to the date of application of the powder and liquid. Do not use after their expiration dates.
- Polisol universal expansion liquid is suitable for rapid or slow heating.
- The investment material contains silica powder. Inhalation of this powder can be harmful and can cause irreversible damage to the lungs (silicosis, lung cancer).
- Store the powder and liquid at a normal room temperature (23°C/73.4°F). If the liquid is exposed to temperatures at or below 5°C (41°F), the liquid will freeze and cannot be used thereafter.
- The investment material should not come into contact with plasters or plaster investment materials.
- The instructions are valid for temperatures of the mass and the room between 23°C (73.4°F)

MODEL CASTING TECHNIQUE

MIXING POWDER AND LIQUID

Mixing ratio: 100 g of powder / 22 ml of liquid
Time of mixing in the vacuum mixer: 60 seconds
Working time: 4 - 5 minutes

DUPLICATION AND INSERTION

Table with 5 columns: MODEL, POWDER, LIQUID, CONCENTRATION OF EXPANSION LIQUID, Expansion liquid/distilled water. Rows for MODEL, CILINDER, and CILINDER with powder and liquid amounts.

RECOMMENDATION: when duplicating the model using duplicating gel with 150 g of powder (1 model) and 33 ml of liquid (16.5 ml of expansion liquid and 16.5 ml of distilled water), we recommend that you:

- accurately measure the prescribed quantity of expansion liquid and water
- mix the investment material in the vacuum mixer for 100 seconds
- the temperature of the powder and liquid should not be lower than 23°C (73.4°F)

Table with 3 columns: MODEL, DUPLICATING GEL, SILICONE. Rows for Model setting time and Hardening the model.

WARNING: Special caution should be taken when removing the duplicated model from the silicone. Use compressed air to help you.

HEATING

Table with 3 columns: After the beginning of mixing, the cylinder is placed in the furnace for a period of, RAPID, SLOW. Rows for Temperature of the preheated furnace, Heating, and The period of maintaining the final temperature.

WARNING: The temperature of the preheated furnace should not exceed 900°C (1652°F) when the cylinder is inserted. After being inserted into the furnace, the cylinder can be heated up to 1050°C (1922°F).

CROWN AND BRIDGEWORK TECHNIQUE

MIXING THE POWDER AND LIQUID

Mixing ratio: 100 g of powder / 22 ml of liquid

Time for mixing in the vacuum mixer:

60 seconds

Working time: 5 minutes


Table with 3 columns: Cylinder size, POWDER, LIQUID. Rows for x 1, x 3, x 6, x 9.

The higher the concentration of expansion liquid, the higher the expansion of the mass. Avoid increasing the density of the consistency, as this will substantially reduce the working time and increase the expansion of the mass!

Table with 4 columns: Precious and semi-precious alloys with a melting point below 1150°C, Precious and semi-precious alloys FOR CERAMICS, Non-precious alloys Ni-Cr, Non-precious alloys Co-Cr. Rows for Cylinder size and Exp. liquid.

Table with 3 columns: After the beginning of mixing, the cylinder is placed in the furnace for a period of, RAPID, SLOW. Rows for Temperature of the preheated furnace, Heating, and The period of maintaining the final temperature.

WARNING: The temperature of the preheated furnace should not exceed 900°C (1652°F) when the cylinder is inserted. After being inserted into the furnace, the cylinder can be heated up to 1000°C (1832°F).

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