© 2024 PALA®

PALA® Mill&Cast Teeth

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(GB) INSTRUCTIONS FOR USE

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© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

© 2024 PALA®

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(GB) INSTRUCTIONS FOR USE

(GB) INSTRUCTIONS FOR USE

(PALA®) Mill&Cast Teeth

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(PALA®) Mill&Cast Teeth

	Spezifikation	Typischer Durchschnittswert	Methode gemäß
Biegebruchfestigkeit	> 65 MPa // > 50 MPa	90 MPa	EN ISO 20795-2 // EN ISO 10477
Biegemodul	> 2000 MPa	2700 MPa	EN ISO 20795-2
Sorption	< 32 µg/mm³	25 µg/mm³	EN ISO 20795-2
Löslichkeit	< 1,6 µg/mm³	0,6 µg/mm³	EN ISO 20795-2
Restmonomer	< 2,2%	0,5%	EN ISO 20795-2 EN ISO 20795-1
Vickers-Härte	> 24 HV	26 HV	/

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Dimensions: full height of disc: 35 mm // tooth part height: 15mm //diameter: 98,5 mm

Shades: A2–A3–A3,5

INTENDED PURPOSE:
PALA® Mill&Cast Teeth discs are suitable for milling full tooth arch to manufacture denture. This medical device is part of the PALA® Mill&Cast Workflow and is always used in combination with Kulzer approved denture base materials to mill dentures.

COMPOSITION: PMMA material with embedded pigments and residual methyl methacrylate less than 1%

INDICATION: Edentulism

CONTRAINDICATION: Do not use the product in case of known allergies or hypersensitivity to acrylic materials.

NOTE! Application on implants is not recommended due to the increased risk of fractures. When implant restoration is required, the screw channel may further increase the risk of fractures.

WEARING TIME: The recommended wearing time is up to 5 years.

INSTRUCTIONS FOR USE
NOTE! Before use of the disc carefully read the label of the disc to choose the correct shade.
Procedure: For proper procedure of design and milling the full tooth arch, positioning the disc and preventing it from moving from the zero point, follow the instruction for use included in the denture base material PalaXpress® manufactured by Kulzer GmbH and PALA® Mill&Cast application guide. The instruction for use PalaXpress® provides also instruction for dosing, handling and polymerization of denture base material. Milling of denture base from the base material are provided in PALA® Mill&Cast application guide.
NOTE! Make sure to thoroughly clean the restoration to remove any residue of the milling and finishing processes using common dental laboratory practice (e.g. steam jet). For better bonding roughen and wet the bonding surfaces of milled teeth with PalaXpress® monomer.
Finishing
Tools for acrylate are used for finishing and final adjustments. For polishing use pumice stone, pastes and soft polishing brushes to achieve a glossy surface. Pay attention not to overheat the material!

INFORMATION FOR DENTIST: In accordance with standard dental medicine practices, dentist must clean medical device before inserting it into the oral cavity. The dentist must provide the patient with instructions for adequate cleaning of the medical device and follow up the patient at regular appointments to assess the condition of the medical device. We recommend cleaning with soft brushes and mild soap for removable medical device and non-abrasive tooth pastes for fixed restorations.

STORAGE: Keep away from direct sunlight!

DISPOSAL: Disposal should be made in accordance with local, state and national legislation.

WARNINGS!
Professional use only! Dental CAD/CAM skills are required to properly process the product. The discs must be milled using a compatible CAD/CAM system. It is important not to overheat the material! Do not use the material after expiry date. Expiry date is printed on the packaging. PMMA dust may develop during the milling process, leading to mechanical irritation of the eyes, skin, and respiratory tract. Make sure to use the extraction system in the workplace. When processing PMMA discs, wear personal protective equipment (protection mask for dust, safety glasses, ...). To ensure traceability and assignment of the product data as well as the LOT numbers for the patient at any time, it is recommended to store the milling blank in the original packaging during two milling operations.

Further safety-relevant information can be found in the material safety data sheet available on request.

In case that any serious incident has occurred in relation to medical device PALA® Mill&Cast Teeth disc it should be reported to the manufacturer (Polident d.o.o.; vigilanca@polident.si) and the competent authority of the Member State, here the user is established.

	Catalogue number		Keep away from sunlight
	Batch code		Medical device
	Use by date		CE marking
	Consult instructions for use		Unique Device Identification System
	Manufacturer		Distributor

	Especificación	Valor medio típico	Método conforme a
Resistencia a la flexión	> 65 MPa // > 50 MPa	90 MPa	EN ISO 20795-2 // EN ISO 10477
Módulo de flexión	> 2000 MPa	2700 MPa	EN ISO 20795-2
Sorción	< 32 µg/mm³	25 µg/mm³	EN ISO 20795-2
Solubilidad	< 1,6 µg/mm³	0,6 µg/mm³	EN ISO 20795-2
Monómero residual	< 2,2%	0,5%	EN ISO 20795-2 EN ISO 20795-1
Dureza Vickers	> 24 HV	26 HV	/

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