PALA® Mill Base

**KULZER**  
MITSUI CHEMICALS GROUP

Kulzer GmbH  
Leipziger Straße 2  
63450 Hanau (Germany)

POLIDENT d.o.o.  
Voljca Draga 42, 5293 Voljca Draga, Slovenija  
<http://www.polident.si>

VER 27.05.2024\_MDR 01

No. 0200316

**CE 1304**

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**TECHNISCHE PARAMETER**

	Specification	Typical average value	Method according to
Flexural strength	> 65 MPa	90 MPa	EN ISO 20795-1
Flexural modulus	> 2000 MPa	2700 MPa	EN ISO 20795-1
Sorption	< 32 µg/mm³	25 µg/mm³	EN ISO 20795-1
Solubility	< 1,6 µg/mm³	0,6 µg/mm³	EN ISO 20795-1
Residual monomer	< 2,2%	0,5 <span> </span> %	EN ISO 20795-1

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**Dimensions:** thickness: 30mm // diameter: 98,5 mm

**Shades:** PINK – PINK LIVE

**INTENDED PURPOSE:** **PALA® Mill Base** disc is suitable for milling full and partial dentures and temporary immediate load implant supported dentures.

**COMPOSITION:** PMMA material with embedded pigments and residual methyl methacrylate less than 1%

**INDICATION:** Edentulism, partial edentulism

**CONTRAINDICATIONS:** Do not use the product in case of known allergies or hypersensitivity to acrylic materials.

**WEARING TIME:** The recommended wearing time is up to 5 years.

**INSTRUCTIONS FOR USE**

**NOTE!** Before use of the disc carefully read the label of the disc to choose the correct height of the disc and shade.

**Design rules to consider within CAD software for the design of denture bases**

Design the denture base with suitable CAD dental software.

Fix the disc in the milling machine holder. Load the CAM -file on the milling machine. Milling strategies for PMMA must be considered according to the milling machine instructions. The denture is separated from the disc with standard grinding tools for PMMA resin after milling.

**Bonding of denture teeth**

Bond the artificial teeth to the denture base with PMMA resin or light curing composite material according to manufacturer’s instructions. The artificial teeth must be roughened and prepared according to manufacturer’s instructions.

The denture base does not require special treatment in case of PMMA bonding resin but usually requires a primer treatment in case of composite bonding material according to manufacturer’s instructions. The denture can be finished with standard grinding and polishing tools for PMMA resin after the bonding of artificial teeth.

**Finishing**

Tools for acrylate are used for finishing and final adjustments. For polishing use grinding paper, polishing rubbers, pumice stones, pastes and soft polishing brushes suitable for acrylic resins to achieve glossy surface.

**CLEANING:** Make sure to thoroughly clean the restoration to remove any residue of the milling and finishing processes using common dental laboratory practice (e.g. steam jet or mild soap/detergent).

**INFORMATION FOR DENTIST:** In accordance with standard dental medicine practices, dentist must clean medical device before inserting it into the oral cavity. The dentist must provide the patient with instructions for adequate cleaning of the medical device and follow up the patient at regular appointments to assess the condition of the medical device. We recommend cleaning with soft brushes and mild soap for removable medical device and non-abrasive tooth pastes for fixed restorations.

**STORAGE:** Keep away from direct sunlight!

**DISPOSAL:** Disposal should be made in accordance with local, state and national legislation.

**WARNINGS:** Professional use only! Dental CAD/CAM skills are required to properly process the product. The discs must be milled using a compatible CAD/CAM system. It is important not to overheat the material!

Do not use the material after expiry date. Expiry date is printed on the packaging.

PMMA dust may develop during the milling process, leading to mechanical irritation of the eyes, skin, and respiratory tract. Make sure to use the extraction system in the workplace. When processing PMMA discs, wear personal protective equipment (protection mask for dust, safety glasses, ...).

To ensure traceability and assignment of the product data as well as the LOT numbers for the patient at any time, it is recommended to store the milling blank in the original packaging during two milling operations.

Further safety-relevant information can be found in the material safety data sheet available on request.

In case that any serious incident has occurred in relation to the device it should be reported to the manufacturer (Polident d.o.o.; [vigilancia@polident.si](mailto:vigilancia@polident.si)) and the competent authority of the Member State, where the user is established.

**List of symbols**

	Catalogue number		Keep away from sunlight
	Batch code		Medical device
	Use by date		CE marking
	Consult instructions for use		Unique Device Identification System
	Manufacturer		Distributor

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