

PoliDent

🏠 *since 1958*



ACRYLICS

ACRYLICS

PoliDent



POLIDENT d.o.o. Dental Products Industry

Volčja Draga 42
Si-5293 Volčja Draga
Slovenia

Phone: +386 5 330 48 40
Fax: +386 5 330 48 70
E-mail: polident@polident.si

www.polident.si

www.polident.si

Date: 24.01.2019

poliHOT

Acrylate for the fabrication of full and partial dentures by using the heat curing

Characterictis:	Colour matching with the acrylate for the fabrication of repairs (Policold).
Indications:	full dentures, partial dentures
Color range:	pink, pink veined, M26, M26 V
Packing:	02866 Polihot powder pink - 700 g 02876 Polihot powder pink veined - 700 g 03064 Polihot powder M26 - 700 g 03065 Polihot powder M26V - 700 g 02867 Polihot liquid - 500 g

*for larger orders we can arrange other quantities of packaging



poliCOLD

Acrylate for the fabrication of partial dentures, repairs and rebasings by using the cast technique and the cold curing

Characterictis:	Colour matching with acrylate for the fabrication of dentures using the heat curing (Polihot).
Indications:	partial denture, repairs, rebasings
Color range:	pink, pink veined
Packing:	02868 Policold powder pink - 700 g 02879 Policold powder pink veined - 700 g 02869 Policold liquid - 500 g 02875 Complete: Policold powder pink - 200 g, Policold liquid - 150 g 02885 Complete: Policold powder pink veined - 200 g, Policold liquid - 150 g

*for larger orders we can arrange other quantities of packaging



K-30 S

Acrylate for the fabrication of full and partial dentures by using the packing / pressing technique and the classical heat curing

Characterictis:	The use of the cold curing material Polirepar S is recommended, the colours of which matches the acrylate for the fabrication of prothesis by using the heat curing (K-30 S).
Indications:	full dentures, partial dentures
Color range:	pink, pink veined, PN, PN veined, M26, M26V
Packing:	02706 powder pink - 700 g 02719 powder pink veind -700g 02721 powder PN - 700 g 02722 powder PN veined - 700 g 02734 powder M26V - 700g 02736 powder M26 - 700g 02707 liquid - 500 g

*for larger orders we can arrange other quantities of packaging



poliREPAR S

Acrylate for the fabrication of repairs by using the cold curing

Characterictis:	Colour matching with acrylate for the fabrication of dentures using the heat curing K-30S.
Indications:	repairs
Color range:	pink, PN, PNV
Packing:	02948 Polirepar S powder - 200 g 03094 Polirepar S powder PN - 200 g 03095 Polirepar S powder PNV - 200 g 02716 Polirepar S liquid - 150 g

*for larger orders we can arrange other quantities of packaging



ORTOpoli

Acrylate for the fabrication of orthodontic appliances by using the spray-on technique and the cold curing.

Characterictis:	Fast mixing of the powder with the liquid, the material therefore doesn't spread all over the model. Large choice of basic colours. The colouring agents are cadmium free, therefore it isn't toxic.
Indications:	orthodontic appliance
Color range:	red, yellow, blue, colourless, fl. yellow, orange, indigo blue, green, pink
Packing:	02872 powder, colourless - 700 g, shaker 02873 liquid, red - 500 g 02882 liquid, yellow - 500 g 02883 liquid, blue - 500 g 02884 liquid, colourless - 500 g 02979 liquid, fl. yellow - 150 g 02980 liquid, orange - 150 g 02981 liquid, indigo blue - 150 g 02982 liquid, green - 150 g 02983 liquid, pink - 150 g

*for larger orders we can arrange other quantities of packaging



poliTRAY

Acrylate for the fabrication of individual impression trays by using the cold curing

Characterictis:	Short polymerisation. Microporous surface enables a good hold to the impression material.
Indications:	Individual impression trays
Color range:	blue
Packing:	02870 Politray powder - 700 g 02871 Politray liquid - 500 g 02997 Politray liquid - 150 g

*for larger orders we can arrange other quantities of packaging



uniTRAY LC

Light-curing material in form of prefabricated, upper baseplates to make rigid, custom impression trays and bite-blocks. This material can be cured in all ultraviolet and Halogen light curing units.

Characterictis:	Place the tray with the model in a light curing unit for 6 - 10 min.
Indications:	custom impression trays, bite-blocks
Color range:	pink
Packing:	15000 Unitray LC - box with 50 upper baseplates.



poliTEMP

Cold – curing acrylic resin for temporary crowns and bridges

Characterictis:	High precision and thus a better contact of the crown with the gingival area of the core. Glossy surface and minor accumulation of plaque because of the non porosity of the final product. Highly abrasion resistant. Simple correction. Not toxic.
Indications:	temporary crowns, temporary bridges
Color range:	A1, A2, A3, A3,5, B3, C2, D3, Incisal
Packing:	02904 powder A1 - 80 g 02949 powder A2 - 80 g 02888 powder A3 - 80 g 02889 powder A 3,5 - 80 g 02890 powder B3 - 80 g 02891 powder C2- 80 g 02892 powder D3 - 80 g 02893 powder incisal - 80 g 02894 liquid - 150 g

