

1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND COMPANY/UNDERTAKING 1.1 Product identifier Trade name CAD-CAM MODEL DISC 1.2 Relevant identifier uses of the substance or mixture and uses advised against Recommended use Milling of stumps/dies, working models and antagonists. Verifying frameworks on plaster model before final ZrO ₂ milling. 1.3 Details of the supplier of the Safety Data Sheet Manufacturer/supplier Polident d.o.o., Dental Products Industry Volčja Draga 42, 5293 Volčja Draga, Slovenija Tel.:00386 5 3304840, Fax: 00386 5 3304870 E-mail: polident@polident.si Information provided by 00386 5 3304859 1.4 Emergency telephone number 00386 5 3304840	
2. HAZARDS IDENTIFICATION 2.1 Classification of the substance or mixture Classification according to Regulation 1272/2008/EC: Not applicable 2.2 Label elements Labelling according to Regulation 1272/2008/EC: Not applicable 2.3 Other hazards None Results of PBT and vPvB assessment PBT: Not applicable. vPvB: Not applicable.	
3. COMPOSITION/INFORMATION ON INGREDIENTS 3.1 Chemical characterization: Mixtures Description Product is composed of polyurethane. Hazardous ingredients according to Regulation 1272/2008/EC: Not applicable	
4. FIRST AID MEASURES 4.1 Description of first aid measures General information No special measures required. Inhalation Remove the casualty into fresh air. In case of complaints get medical advice. Skin contact Wash off with water and soap. Consult a physician if irritation persists. Eye contact Rinse eye with plenty of water if it is necessary seek medical attention: only mechanical effect.	

	Ingestion Wash out mouth with water and drink water. In case of complaints get medical advice.
4.2	Most important symptoms and effects, both acute and delayed None known.
4.3	Indication of any immediate medical attention and special treatment needed Treat symptomatically.
5.	FIRE – FIGHTING MEASURES
5.1	Extinguishing media Suitable extinguishing media Water, foam, dry powder, carbon dioxide. Unsuitable extinguishing media for safety reasons None.
5.2	Special hazards arising from the substance or mixture Dangerous gases may be released in case of fire.
5.3	Advice for firefighters Wear self-contained breathing apparatus.
6.	ACCIDENTAL RELEASE MEASURES
6.1	Personal precautions, protective equipment and emergency procedures Not required.
6.2	Environmental precautions No special measures required.
6.3	Methods and material for containment and cleaning up Pick up mechanically. Dispose of in accordance with regulations.
6.4	Reference to other sections For information on safe handling see Section 7. For information on personal protection equipment see Section 8. For information on disposal information see section 13.
7.	HANDLING AND STORAGE
7.1	Precautions for safe handling Safe handling advice Only adequate trained persons may deal with product. For use in dentistry only. Keep away from children. Ensure the area is well ventilated. No eating, drinking or smoking during handling with product. Follow the usual good standards of occupational hygiene. Information on fire and explosion protection No special measures required.
7.2	Conditions for safe storage, including any incompatibilities Storage Requirements for storage areas Keep in dry storage.
7.3	Specific end use(s) None.
8.	EXPOSURE CONTROLS /PERSONAL PROTECTION
8.1	Control parameters Ingredients with limit values that require monitoring at the workplace The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
8.2	Exposure controls General protective measures Provide adequate ventilation. Hygiene measures Follow the usual good standards of occupational hygiene and dental practise.

Respiratory protection

Avoid contact with powder generated by final working.

Wear respiratory protective device in case of dust formation.

Hand protection

Not required.

Eye protection

Wear safety goggles during mechanical processing.

9. PHYSICAL AND CHEMICAL PROPERTIES

9.1 Information on basic physical and chemical properties

Appearance:

Form: firm disc

Colour: brown and beige

Odour: odourless

Changes in condition: /

Melting point: Undetermined

Boiling point: Undetermined

Flash point: >250°C

Ignition temperature: >350°C

Lower explosion limit: 1,1 vol%

Upper explosion limit: 8 vol%

Vapour pressure: not applicable

Density: 1,0 g/cm³ (beige), 0,7 g/cm³ (brown)

Viscosity (dynamic): not applicable

Solubility in water: insoluble

pH – value: not applicable

9.2 Other information

None.

10. STABILITY AND REACTIVITY

10.1 Reactivity

No relevant information available.

10.2 Chemical stability

The product is stable under normal handling and storage conditions.

10.3 Possibility of hazardous reactions

No hazardous reactions known.

10.4 Conditions to avoid

Decompositions begins at appr. 250°C.

10.5 Incompatible materials

No relevant information available.

10.6 Hazardous decomposition products

None when used as prescribed.

Dangerous vapours may be released in case of fire: carbon monoxide, carbon dioxide, nitrogen oxide, hydrogencyanide, etc.

11. TOXICOLOGICAL INFORMATION

11.1 Information on toxicological effects

Acute oral toxicity

No toxiciting effects known.

Irritation

No irritating effects known.

Sensitization

No sensitizing effects known.

12. ECOLOGICAL INFORMATION
12.1 Toxicity Hazardous to the aquatic environment: No further relevant information available, no evidence for hazardous properties
12.2 Persistence and degradability: No further relevant information available, no evidence for hazardous properties
12.3 Bioaccumulative potential: No further relevant information available, no evidence for hazardous properties
12.4 Mobility in soil: No further relevant information available, no evidence for hazardous properties
12.5 Results of PBT and vPvB assessment PBT: Not applicable. vPvB: Not applicable.
12.6 Other adverse effects: No further relevant information available.
13. DISPOSAL CONSIDERATIONS
13.1 Waste treatment methods Product: Disposal should be made in accordance with local, state and national legislation. Uncleaned packaging: Disposal should be made in accordance with local, state and national legislation. Code of waste EWC: 18 01 07
14. TRANSPORT INFORMATION
14.1 UN Number: see section 14.2
14.2 UN proper shipping name Land transport ADR: Not dangerous according to transport regulations. Land transport RID: Not dangerous according to transport regulations. Inland waterway transport ADN: Not dangerous according to transport regulations. Shipment by sea IMDG: Not dangerous according to transport regulations. Air transport ICAO/IATA: Not dangerous according to transport regulations.
14.3 Transport hazard class(es): see section 14.2
14.4 Packing group: see section 14.2
14.5 Environmental hazards: None.
14.6 Special precautions for user: see section 14.2
14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code: None.
14.8 Transport additional information: None.
15. REGULATORY INFORMATION
15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture National regulations Other regulations, limitations and prohibitive regulations
16. OTHER INFORMATION This information is based on our present knowledge and experience. They are intended as a description of the products safety requirements and are not to be seen as a quarantee of certain product`s features. New versions of Safety Data Sheet will be <i>laying printed</i> for the data compared to the previous version altered.